

Hu. Coagulation Factor XI ELISA

Product Number: CF46

Store at 4°C

FOR RESEARCH USE ONLY

Document Control Number: CF46.110629

Page 1 of 4

Total Human Coagulation Factor XI ELISA

For Research Use Only

INTRODUCTION

This assay is intended for the quantitative determination of total Factor XI antigen in human plasma.

Factor XI is a disulfide linked two-chain glycoprotein zymogen and is the precursor of the coagulation enzyme Factor XIa¹. Factor XI consists of two identical monomers and circulates in plasma in complex with kininogen². Factor XI is activated by Factor XIIa and converts Factor IX to Factor IXa during the intrinsic pathway of the coagulation cascade³.

PRINCIPLES OF PROCEDURE

Human Factor XI will bind to the affinity purified capture antibody coated on the microtiter plate. Factor XI and XIa will react with the antibody on the plate. After appropriate washing steps, peroxidase labeled polyclonal anti-human Factor XI primary antibody binds to the captured protein. Excess antibody is washed away and TMB substrate is used for color development at 450nm. A standard calibration curve is prepared along with the samples to be measured using dilutions of human Factor XI. Color development is proportional to the concentration of Factor XI in the samples.

MATERIALS PROVIDED

Component	Contents	Quantity	Storage	Cat. No.
Coated Plate	Anti-human Factor XI coated 96-well plate	1 plate	4°C	CF46a
Standard	Lyophilized Human Factor XI standard	1 vial	4°C	CF46b
Wash Buffer	10x solution for washing plate	50 mL	4°C	CF46c
Primary Antibody	Lyophilized Anti-human Factor XI-HRP Conjugate	1 vial	4°C	CF46d
Substrate	TMB Substrate	10 mL	4°C	CF46e

MATERIALS NEEDED BUT NOT PROVIDED

- 1. Pipettes covering 0-10 μ l and 200-1000 μ l and tips
- 2. 12-channel pipette covering 30-300 μ l
- 3. 1 N H₂SO₄
- 4. DI water
- 5. Microtiter plate spectrophotometer with a 450 nm filter
- 6. Microtiter plate shaker with uniform horizontally circular movement up to 300 rpm

STORAGE CONDITIONS

- 1. Store this kit and its components at 4°C until use.
- 2. The reconstituted Standard and Primary Antibody may be stored at -70°C for later use. **DO NOT** freeze/thaw the Standard or Primary Antibody more than once.

PROCEDURAL NOTES

- 1. Use aseptic technique when opening and dispensing reagents.
- 2. This kit is designed to work properly as provided and instructed. Additions, deletions, or substitutions to the procedure or reagents are not recommended, as they may be detrimental to the assay.
- 3. Exercise universal precautions during the performance or handling of this kit or any component contained therein.

SAMPLE COLLECTION AND PREPARATION

The assay measures total human Factor XI in the 0.1-100 ng/ml range. Samples giving human Factor XI levels above 100 ng/ml should be diluted in Blocking Buffer before use. A 1:1,000 to 1:5,000 dilution for plasma is suggested for best results.

REAGENT PREPARATION

- 1. **10x Wash Buffer:** Dilute the 50 mL of concentrate to 1x with 450 mL of DI water prior to use.
- 2. **TBS Buffer:** 0.1 M Tris, 0.15 M NaCl, pH 7.4.
- 3. **3% BSA Blocking Buffer:** 3% BSA in TBS Buffer.
- 4. **Standard:** Reconstitute with DI Water as directed on the vial to give a 1000 ng/mL Standard Stock Solution. Prepare immediately prior to use.
- 5. **Primary Antibody:** Reconstitute with 3% BSA Blocking Buffer as directed on the vial and vortex gently to mix. Prepare immediately prior to use.

SANDARD PREPARATION

Prepare the Standard Stock Solution as described above. Do not prepare the standards until you are ready to apply them to the plate.

Table 1: Preparation of Standard Curve

Standard	Factor XI Concentration (ng/mL)	Blocking Buffer (µL)	Transfer Volume (µL)	Transfer Source	Final Volume (µL)	
S9	100	900	100	Stock Vial	500	
S8	50	500	500	S9	600	
S ₇	20	600	400	S ₈	500	
S ₆	10	500	500	S ₇	500	
S ₅	5	500	500	S ₆	600	
S4	2	600	400	S ₅	500	
S3	1	500	500	S ₄	500	
S ₂	0.5	500	500	S ₃	600	
S ₁	0.2	600	400	S ₂	1000	
S ₀	0	500			500	

ASSAY PROCEDURE

- 1. Add 100 μ l of the Standards and unknowns to the wells in duplicate. Shake the plate at 300 rpm for 30 minutes at room temperature (RT). See Scheme I for a suggested plate layout.
- 2. Wash the plate 3 times according to the following wash procedure:
 - a. Remove the contents of each well by inversion of the plate.
 - b. Tap out the remaining contents of the plate onto a lint free paper towel.
 - c. Add 300 μ L of 1x Wash Buffer.
 - d. Let stand for 2-3 minutes.
 - e. Repeat procedure two more times, then proceed to step "f".
 - f. Remove the contents of each well by inversion of plate into an appropriate disposal device.
 - g. Tap out the remaining contents of the plate onto a lint free paper towel, then proceed to step 3.
- 3. Add 100 µl of the Primary Antibody to each well. Shake the plate at 300 rpm for 30 minutes at RT.
- 4. Wash the plate three times as in step 2.
- 5. Add 100 µl of TMB Substrate to each well. Shake the plate at 300 rpm for 5-10 minutes at RT.
- 6. Stop the reaction by adding 50 μ l of 1N H_2SO_4 to each well and read the plate at 450 nm.


Scheme .	I: S	uggesi	ted	Plate	Layo	out
_		_	_			_

	1	2	3	4	5	6	7	8	9	10	11	12
A	S9	S ₈	S7	S ₆	S ₅	S4	S3	S ₂	S ₁	S ₀	U_1	U ₂
В	S9	S ₈	S 7	S ₆	S ₅	S ₄	S 3	S ₂	s_1	S_0	U_1	U_2
C	U3	U ₄	U ₅	U ₆	U7	U ₈	U9	U ₁₀	U ₁₁	U ₁₂	U13	U14
D	U3	U4	U_5	U_6	U7	U_8	U9	U_{10}	U_{11}	U_{12}	U13	U14
		U16										
F	U_{15}	U_{16}	U_{17}	U_{18}	U19	U_{20}	U_{21}	U_{22}	U_{23}	U_{24}	U_{25}	U_{26}
		U_{28}										
Н	U27	U_{28}	U29	U30	U31	U32	U33	U34	U35	U36	U37	U38

CALCULATIONS

- 1. Plot the A_{450} against the concentration of Factor XI in the standards.
- 2. Fit a straight line through the points using a linear fit procedure.
- 3. Calculate the Factor XI concentrations in the unknowns using the equation generated by the standard curve.

Figure 1: Typical Standard Curve Human Coagulation Factor XI



EXPECTED VALUES

The concentration of Factor XI in normal human plasma ranges from 3.0 to 6.0 μg/ml⁴.

REFERENCES

- 1. Kurachi, K, and Davie, E.W.; (1981) Methods in Enzymology: 80: 211-220
- 2. Thompson, R.E., et al.; (1977) J. Clin, Invest.: **60**: 1376
- 3. Walsh, P.N.; (1993) Methods in Enzymology: **222**:65-96
- 4. Bouma, B.N., et al.; (1983) Blood: **62**:1123-1131

DISCLAIMER

This information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. Oxford Biomedical Research, Inc. shall not be held liable for any damage resulting from handling or from contact with the above product. See catalog for additional terms and conditions of sale.

TECHNICAL SUPPORT

If you need technical information or assistance with assay procedures, please call our Technical Support Department at 800-692-4633 or 248-852-8815. Our staff will be happy to answer your questions about this or any other product in the Oxford Biomedical line.

GUARANTEE AND LIMITATION OF REMEDY

Oxford Biomedical Research, Inc. makes no guarantee of any kind, expressed or implied, which extends beyond the description of the material in this ELISA kit, except that these materials and this kit will meet our specifications at the time of delivery. Buyer's remedy and Oxford Biomedical Research, Inc.'s sole liability hereunder is limited to, at Oxford Biomedical Research, Inc.'s option, refund of the purchase price of, or the replacement of, material that does not meet our specification. By acceptance of our products, Buyer indemnifies and holds Oxford Biomedical Research, Inc. harmless against, assumes all liability for the consequence of its use or misuse by the Buyer, its employees, or others. Said refund or replacement is conditioned of Buyer notifying Oxford Biomedical Research, Inc. within thirty (30) days of the receipt of product. Failure of Buyer to give said notice within said thirty (30) days shall constitute a waiver by the Buyer of all claims hereunder with respect to said material(s).

Oxford Biomedical Research, Inc. P.O. Box 522 Oxford, MI 48371 U.S.A.

Orders: 800-692-4633 Technical Service: 248-852-8815 Fax: 248-852-4466 E-mail: info@oxfordbiomed.com

Made in the U.S.A.