

Product Specification Sheet

Description: CypExpress™ 2B6
Product Number: CE2B6.10
Aliquot: 10.0 Gram
Lot Number: ce2B6.200120
Storage: -70°C

DESCRIPTION:	CypExpress™ 2B6 is a permeabilized and stabilized dried yeast powder preparation containing recombinant human CYP 2B6 and recombinant human P450 NADPH oxidoreductase.
ACTIVITY:	A reaction mixture containing 20 µM 7-ethoxy-4-(trifluoromethyl)coumarin and 100 mg/mL of CypExpress™ 2B6 produced 2.6 uM 7-Hydroxy-4-(trifluoromethyl)coumarin in 3 hours at 30°C under the conditions described.
THAWING:	It is important to avoid water adsorption during thawing. Remove the CypExpress™ vial from storage and warm to room temperature prior to opening. Return unused powder to -70°C.
STORAGE:	Store at -70°C. Although CypExpress™ preparations can be warmed to room temperature multiple times, it is best to aliquot as needed to avoid multiple freeze/thaw cycles.
SUGGESTED PROTOCOLS:	<u>The basic protocols attached are intended as a starting point that can be optimized by an experienced user for specific applications and samples.</u>