



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Prothrombin (Mouse) IgG Fraction

Product Number: CT04

Aliquot: 1.0 mg

Lot Number: Typical Lot

Storage: -70°C

DESCRIPTION:	This polyclonal sheep antibody is specific for mouse prothrombin/thrombin and is provided as a purified IgG fraction.
MOLECULAR WEIGHT:	N/A
CONCENTRATION:	10 mg/mL
SIZE:	1.0 mg
VOLUME:	0.1 mL
STORAGE:	-70°C – Avoid repeated freeze/thaw cycles.
FORM:	Frozen liquid
SOURCE:	Sheep serum
BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
PURIFICATION:	IgG fraction purified by immobilized Protein G
SPECIES SPECIFICITY:	Mouse

For research purposes only.