

**Product Specification Sheet**

Anti-Mouse Fibrinogen Biotin  
Labeled IgG  
Product Number: CT 23  
Aliquot: 1.0 mg  
Lot Number: ct23.408  
Storage: -80°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is specific for mouse fibrinogen.
<b>MOLECULAR WEIGHT:</b>	330 kDa
<b>CONCENTRATION:</b>	9.7 mg/mL determined spectrophotometrically.
<b>TAG:</b>	8.0 mole biotin/ mole IgG
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; 1 mM EDTA; pH 6.6
<b>STORAGE:</b>	-80°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Purified mouse fibrinogen from plasma.
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Purified by Protein A Sepharose.
<b>PURITY:</b>	≥ 98% assessed by visual inspection of a Coomassie <sup>®</sup> blue-stained SDS-PAGE gel.
<b>APPLICATIONS:</b>	Suitable for use in Western blot (10 µg/mL) and ELISA (1 µg/mL).
<b>SPECIES SPECIFICITY:</b>	Mouse