

Product Specification Sheet

Description: Bovine Fibrinogen,
FITC Labeled
Product Number: CT51
Aliquot: 10 mg
Lot Number: ct51.1110
Storage: -70°C

DESCRIPTION:	FITC labeled bovine fibrinogen. Prepared from fresh bovine plasma. Plasminogen has been depleted by lysine affinity chromatography. Thaw vial at 37°C without shaking. Mix gently before using.
MOLECULAR WEIGHT:	330 kDa
CONCENTRATION:	13.2 mg/mL
STORAGE BUFFER:	0.02M Sodium citrate-HCL, pH 7.4
STORAGE:	-70°C, Protect from light.
PURITY:	> 95% by SDS-PAGE analysis