



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Sheep Anti-Human tPA Antiserum

Product Number: PA36

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

<b>DESCRIPTION:</b>	Anti Human tPA Polyclonal antiserum (host sheep)
<b>MOLECULAR WEIGHT:</b>	N/A
<b>CONCENTRATION:</b>	Whole antiserum
<b>SIZE:</b>	1.0 mL
<b>VOLUME:</b>	1.0 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum

For research purposes only.