



OXFORD BIOMEDICAL RESEARCH

### Product Specification Sheet

Sheep Anti-Mouse & Rat tPA IgG Fraction  
Product Number: PA39  
Aliquot: 1.0 mg  
Lot Number: Typical Lot  
Storage: -70°C

<b>DESCRIPTION:</b>	Anti-tPA (Mouse & Rat) Sheep Polyclonal, IgG Fraction
<b>MOLECULAR WEIGHT:</b>	160,000
<b>CONCENTRATION:</b>	10 mg/mL
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	0.1 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum
<b>BUFFER:</b>	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>PURIFICATION:</b>	IgG fraction purified by immobilized Protein G.

For research purposes only.