



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Sheep Anti-Rabbit Antiserum

Product Number: PA40

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

<b>DESCRIPTION:</b>	Anti-tPA (Rabbit), Sheep Polyclonal Antiserum
<b>MOLECULAR WEIGHT:</b>	160,000
<b>CONCENTRATION:</b>	Whole antiserum
<b>SIZE:</b>	1.0 mg
<b>VOLUME:</b>	1.0 mL
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Sheep serum
<b>BUFFER:</b>	N/A

For research purposes only.