



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Non-Enzymatic Non-Cleavable Human tPA
Product Number: PA46
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	Human tPA with a double mutation that is permanently single chain and lacks catalytic activity. Recombinant, produced in insect cells.
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	2.0 mg/mL
SIZE:	0.1 mg
VOLUME:	0.05 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture
BUFFER:	0.4M HEPES; 0.1M NaCl; pH 7.4

For research purposes only.