



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Mouse tPA, Non-Enzymatic
Product Number: PA76
Aliquot: 0.1 mg
Lot Number: Typical Lot
Storage: -70°C

DESCRIPTION:	This single chain mouse tPA has the active site serine mutated to alanine rendering the enzyme catalytically inactive. The tPA still retains exosite binding as well as other biological properties of tPA including but not limited to surface binding to fibrinogen.
MOLECULAR WEIGHT:	70,000
CONCENTRATION:	2.0 mg/mL
SIZE:	0.1 mg
VOLUME:	0.05 mL
STORAGE:	-70°C
FORM:	Frozen liquid
SOURCE:	Insect cell culture
BUFFER:	0.4M HEPES; 0.1M NaCl; pH 7.4

For research purposes only.