



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Sheep Anti-Human PAI-1 Antiserum

Product Number: PI57

Aliquot: 1.0 mL

Lot Number: Typical Lot

Storage: -70°C

DESCRIPTION: Anti-PAI-1 (human), sheep polyclonal antiserum.

MOLECULAR WEIGHT: N/A

CONCENTRATION: Whole antiserum

SIZE: 1.0 mL

VOLUME: 1.0 mL

STORAGE: -70°C

FORM: Frozen liquid

SOURCE: Sheep serum

BUFFER: N/A

For research purposes only.