



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Human PAI-1 MoAB

Product Number: PT 10

Aliquot: 1.0 mg

Lot Number pt10.205A

Storage: -80°C

DESCRIPTION:	This is a mouse monoclonal antibody is specific for human plasminogen activator inhibitor-1 (PAI-1) in its complexed form with tissue plasminogen activator (tPA) and urokinase (uPA).
MOLECULAR WEIGHT:	50 kDa glycosylated 43 kDa non-glycosylated
CLONE:	3E529
CONCENTRATION:	4.4 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.01 M Sodium Phosphate; 1 mM EDTA; pH 6.6
STORAGE:	-80°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Complexed PAI-1/tPA
CLASS:	IgG
PURIFICATION:	Purified by Protein A Sepharose.
APPLICATIONS:	This antibody is suitable for use in Western blotting and ELISA.
SPECIES SPECIFICITY:	Human