

Product Specification Sheet

Mouse Fibrinogen
Product Number: CT15
Aliquot: 1.0 mg
Lot Number ct15.316
Storage: -70°C

DESCRIPTION:	This is fibrinogen purified from plasma using a multi-step precipitation process followed by SPE extraction with immunoaffinity resins.
MOLECULAR WEIGHT:	~340 kDa
CONCENTRATION:	1.74 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.02 M Sodium Citrate-HCl; pH 7.4
STORAGE:	-70°C AVOID REPEATED FREEZE THAW CYCLES
PURIFICATION:	Purified using a multi-step precipitation process followed by SPE extraction with immunoaffinity resins.
PURITY:	>95% by SDS-PAGE analysis.
NOTES:	This product is 91% clottable and is homogenous on 4-20% SDS-PAGE gels.