

Product Specification Sheet

Description: 12-LOX Platelet
Standard
Product Number: LX12PS
Aliquot: 200 uL
Lot Number: 063094SS
Storage: -70°C

DESCRIPTION:	12-Lipoxygenase Platelet Standard
MOLECULAR WEIGHT:	~75 kDa
CONCENTRATION:	1.8 mg/mL
STORAGE BUFFER:	25mM Tris, 25 mM potassium phosphate, 300 uM DDC, 20% glycerol, pH 7.8
STORAGE:	Store at -70°C. Avoid repeated freeze/thaw cycles.
SOURCE:	Human platelets.
APPLICATIONS:	Designed for use as a Western blotting standard. Mix 2:1 with SDS sample buffer. Loading 30 uL of the 2:1 mixture is recommended for mini gels.