

Product Specification Sheet

Mouse Anti-Human GST M1-1
Product Number: GS17
Aliquot: 100 µg
Lot Number: gs17.150415
Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody to human glutathione S-transferase M1-1 is protein A & A-G affinity purified and has been tested for cross reactivity against all 18 cytosolic GSTs showing only slight cross reactivity with GSTs M4-4 and M5-5.
CONCENTRATION:	1.00 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN ₃
STORAGE:	-70°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Recombinant human GST M1-1 expressed in <i>E. coli</i>
PURIFICATION:	Protein A & A-G affinity purified
APPLICATIONS:	Enzyme Immunoassay & Immunohistochemistry. Not applicable for Western Blotting (SDS or Native).
SPECIES SPECIFICITY:	Human. Slight cross reactivity with GST M4-4 and M5-5.