

Product Specification Sheet

Mouse Anti-Human GST T1-1
Product Number: GS20
Aliquot: 100 µg
Lot Number: gs20.160509
Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody to human glutathione S-transferase T1-1 is protein A & A-G affinity purified and shows no cross reactivity for all 18 cytosolic GSTs.
CONCENTRATION:	1.00 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN ₃
STORAGE:	-70°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Recombinant human GST T1-1 expressed in <i>E. coli</i>
PURIFICATION:	Protein A & A-G affinity purified
APPLICATIONS:	Enzyme Immunoassay & Immunohistochemistry. Not applicable for Western Blotting (SDS or Native).
SPECIES SPECIFICITY:	Human