

Product Specification Sheet

Mouse Anti-Human GST K1-1
Product Number: GS25
Aliquot: 100 µg
Lot Number: gs25.150209
Storage: -70°C

DESCRIPTION:	This mouse monoclonal antibody to human glutathione S-transferase K1-1 is protein A & A-G affinity purified and shows no cross reactivity for all 18 cytosolic GSTs.
CONCENTRATION:	1.00 mg/mL
STORAGE BUFFER:	Phosphate Buffered Saline, pH 7.4 and 25% Glycerol and 0.02% NaN ₃
STORAGE:	-70°C. Avoid repeated freeze thaw cycles
IMMUNOGEN:	Recombinant human GST K1-1 expressed in <i>E. coli</i>
PURIFICATION:	Protein A & A-G affinity purified
APPLICATIONS:	Enzyme Immunoassay. Immunohistochemistry. Western Blotting: 1:4K dilution is reactive with ≥ 0.2µg of recombinant human GST K1-1 (OBR product #GS71) in conjunction with an HRP conjugated anti-mouse secondary antibody.
SPECIES SPECIFICITY:	Human

