

**Product Specification Sheet**

Anti-Human tPA Antisera  
Product Number: PA57  
Aliquot: 1.0 mL  
Lot Number pa57.406  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit, is provided as neat antiserum, and is specific for human tPA.
<b>MOLECULAR WEIGHT:</b>	70 kDa
<b>CONCENTRATION:</b>	Whole Antiserum
<b>STORAGE BUFFER:</b>	Not Applicable
<b>STORAGE:</b>	-70°C. AVOID REPEATED FREEZE THAWS
<b>IMMUNOGEN:</b>	Recombinant full length wt-Human tPA (Glycosylated).
<b>PURIFICATION:</b>	Not Purified
<b>PURITY:</b>	Not Purified
<b>APPLICATIONS:</b>	Western Blotting: 1:10,000 ELISA: 1:100,000 This antibody is suitable for use in IHC applications.
<b>SPECIES SPECIFICITY:</b>	Human