

**Product Specification Sheet**

Anti-Human uPA IgG Fraction  
Product Number: PA66  
Aliquot: 1.0 mg  
Lot Number pa66.706A  
Storage: -70°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is specific for human uPA. The antibody reacts with both free and complex forms of uPA.
<b>MOLECULAR WEIGHT:</b>	54 kDa High Molecular Weight form 31 kDa Low Molecular Weight form
<b>CONCENTRATION:</b>	14.4 mg/mL determined Spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.15 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Recombinant Human uPA
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Purified by protein A sepharose.
<b>APPLICATIONS:</b>	Suitable for use in ELISA and Western blot.
<b>SPECIES SPECIFICITY:</b>	Human