

Murine PAI-1 (wild type, active)
Product Number: PI16
Aliquot: 0.5 mg
Lot Number: pi16.1117
Storage: -70°C

DESCRIPTION:	This recombinant mouse plasminogen activator inhibitor 1 is a wild type active fraction expressed in <i>E. coli</i> . It typically contains 30-70% active PAI-1.
MOLECULAR WEIGHT:	43 KDa
CONCENTRATION:	1.78 mg/mL
STORAGE BUFFER:	0.05 M Sodium Acetate; 0.10 M NaCl; 1 mM EDTA; pH 5.0
STORAGE:	Store at -70°C and avoid freeze/thaw cycles.
PURITY:	>95% by SDS-PAGE analysis
SOLUBILITY:	> 2 mg/mL and < 5 mg/mL