



Product Specification Sheet

Description: Rat PAI-1
Product Number: PI19
Aliquot: 0.5 mg
Lot Number: PI19.810B
Storage: -70°C

DESCRIPTION: Wild-type active Rat PAI-1 produced in *E. coli*

MOLECULAR WEIGHT: 43 kDa

CONCENTRATION: 4.0 mg/mL

STORAGE BUFFER: 0.05M Sodium Acetate; 0.4M NaCl; 1mM EDTA; pH 5.0

STORAGE: -70°C

PURITY: >95% by SDS-PAGE analysis