

**Product Specification Sheet**

Anti-Human PAI-1 MoAB (non-inhibitory)  
Product Number: PI71  
Aliquot: 1.0 mg  
Lot Number: pi71.605A  
Storage: -70°C

<b>DESCRIPTION:</b>	This monoclonal antibody produced in mouse and specific for human PAI-1 and exhibits non-inhibitory properties. This antibody will react will both free and complexed forms of PAI-1
<b>MOLECULAR WEIGHT:</b>	43 kDa non-glycosylated PAI-1 50 kDa glycosylateed PAI-1
<b>CLONE:</b>	3A120
<b>CONCENTRATION:</b>	5.52 mg/mL determined spectrophotometrically.
<b>STORAGE BUFFER:</b>	0.05 M Sodium Phosphate; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Recombinant Human PAI-1
<b>CLASS:</b>	IgG
<b>PURIFICATION:</b>	Purified by Protein A Sepharose.
<b>APPLICATIONS:</b>	Suitable for use in Western blotting.
<b>SPECIES SPECIFICITY:</b>	Human.