

**Product Specification Sheet**

Description: Chicken Plasminogen  
Product Number: PL40  
Aliquot: 1.0 mg  
Lot Number: pl40.510  
Storage: -70°C

<b>DESCRIPTION:</b>	Chicken plasminogen prepared from fresh plasma.
<b>MOLECULAR WEIGHT:</b>	92 kDa
<b>CONCENTRATION:</b>	2.8 mg/mL
<b>STORAGE BUFFER:</b>	0.1 M HEPES; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C
<b>PURIFICATION:</b>	Purified by immobilized lysine chromatography.
<b>PURITY:</b>	>95% by SDS-PAGE analysis