

**Product Specification Sheet**

Anti-Mouse Antiplasmin MoAB IgG  
(Inhibitory)  
Product Number: PL56  
Aliquot: 1.0 mg  
Lot Number: pl56.503A  
Storage: -70°C

<b>DESCRIPTION:</b>	This mouse monoclonal antibody produced in an antiplasmin knockout mouse. The antibody is provided as an IgG prep and is inhibitory.
<b>MOLECULAR WEIGHT:</b>	67 kDa for mature antiplasmin**
<b>CONCENTRATION:</b>	1.16 mg/mL determined spectrophotometrically.
<b>CLONE:</b>	MAP4H9
<b>STORAGE BUFFER:</b>	0.5 M Sodium Phosphate; 0.1 M NaCl; pH 7.4
<b>STORAGE:</b>	-70°C. <b>AVOID REPEATED FREEZE THAWS.</b>
<b>IMMUNOGEN:</b>	Mouse Antiplasmin
<b>ISOTYPE:</b>	IgG
<b>PURIFICATION:</b>	Purified by Protein A Sepharose.
<b>APPLICATIONS:</b>	Western blotting
<b>SPECIES SPECIFICITY:</b>	Mouse
<b>**REFERENCE:</b>	<i>Hanss, Michel ML a et al.</i> , Blood Coagulation & Fibrinolysis. 14(1):107-111, January 2003.