

Product Specification Sheet

Human Plasminogen, FITC
Labeled
Product Number: PL71
Aliquot: 1.0 mg
Lot Number: p171.919
Storage: -70°C

DESCRIPTION: FITC Labeled Human Plasminogen

MOLECULAR WEIGHT: 92,000

CONCENTRATION: 1.36 mg/mL

STORAGE BUFFER: 0.1 M HEPES, 0.1M NaCl, pH 7.4

STORAGE: -70°C Protect from light.

SOURCE: Human plasma.

LABELING: 2.77 moles FITC per mole protein

PURITY: >95% by SDS-PAGE analysis.