

Product Specification Sheet

Description: Monoclonal Anti-Rat
Cytochrome
Product Number: PM25
Aliquot: 100 µL
Lot Number: 061802DM
Storage: -70°C

| | |
|-----------------------------|--|
| DESCRIPTION: | Monoclonal Anti-Rat Cytochrome P450 (2B1) |
| CONTENTS: | 100 µL per vial. Contains sufficient purified monoclonal mouse IgG for at least 10 Western blots. |
| STORAGE: | Stable for several weeks at 4°C, several months at -20°C, and indefinitely at -70°C. Do not freeze and thaw repeatedly. Refreezing the product may result in diminishment or loss of activity. |
| WORKING DILUTION: | A dilution of 1/1000 for Western blots will yield a strong band with 10 µg of rat microsomes using Alkaline Phosphatase. |
| SPECIES SPECIFICITY: | Identification of 2B1 by immunoassay or Western blots. |