

Product Specification Sheet

Anti-Mouse tPA Antiserum
Product Number: PT30
Aliquot: 1.0 mL
Lot Number pt30.503A
Storage: -70°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is presented as antiserum.
MOLECULAR WEIGHT:	70 kDa
CONCENTRATION:	Not determined
STORAGE BUFFER:	0.05 M Sodium Phosphate; 0.01 M NaCl; 1 mM EDTA; pH 6.6
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Recombinant Mouse tPA
CLASS:	NA
PURIFICATION:	NA
APPLICATIONS:	This antibody exhibits utility in Western applications
SPECIES SPECIFICITY:	Mouse.