

**Product Specification
Sheet**

Receptor Associated Protein
Product Number: RA01
Aliquot: 500 µg
Lot Number ra01.717
Storage: -70°C

DESCRIPTION:	This human receptor associated protein (RAP) was made in <i>E. coli</i> . This RAP binds with high affinity to certain low-density lipoprotein (LDL) receptors such as LRP, gp330/megalin and the very low-density lipoprotein (VLDL).
MOLECULAR WEIGHT:	39 kDa
CONCENTRATION:	1.06 mg/mL determined spectrophotometrically.
STORAGE BUFFER:	0.1M Tris-HCL; 0.15M NaCl; pH 7.4
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
PURITY:	>90% pure by visual inspection of Coomassie™ Blue stained SDS-PAGE gel.
CONSIDERATIONS:	This protein contains endotoxin that should be removed prior to use in active animal or cellular studies.