

**Product Specification Sheet**

Active Recombinant Mouse  
HMW tc-uPA  
Product Number: UP34  
Aliquot: 0.05 mg  
Lot Number: up34.813a  
Storage: -70°C

<b>DESCRIPTION:</b>	Active two-chain HMW mouse urokinase recombinant protein produced in insect cells.
<b>SIZE:</b>	0.05 mg
<b>CONCENTRATION:</b>	1.0 mg/mL
<b>VOLUME:</b>	0.05 mL
<b>MOLECULAR WEIGHT:</b>	44,000
<b>SPECIFIC ACTIVITY:</b>	57,194 IU/mg
<b>BUFFER:</b>	0.05M Sodium Acetate; 0.1M NaCl; 1mM EDTA pH 5.0
<b>STORAGE:</b>	-70°C
<b>FORM:</b>	Frozen liquid
<b>SOURCE:</b>	Insect cell culture
<b>PURITY:</b>	> 95% by SDS-PAGE analysis

For Research Purposes Only.