

Product Specification Sheet

Description: GDF-15 (H6D Variant) Recombinant Human
Product Number: CK01
Aliquot: 20 ug
Lot Number: Typical Lot
Storage: -20°C

DESCRIPTION: Recombinant Human GDF-15 (H6D Variant)

SYNONYMS: MIC-1, Placental TGFbeta, Prostate Differentiation Factor

STORAGE: -20°C

RECONSTITUTION: The lyophilized protein should be reconstituted and aliquoted and stored at -20°C.

SOURCE: Human cell culture

PURITY: > 98% by SDS-PAGE analysis

APPLICATIONS: Cell inhibition assay with DU-145 cells yields an ED₅₀ at 2.0 ug/mL