

Product Specification Sheet

Description: BMP-7,
Recombinant Human Protein
Product Number: CK05
Aliquot: 10 ug
Lot Number: Typical Lot
Storage: -20°C

DESCRIPTION:	Recombinant Human Bone Morphogenetic Protein 7 (BMP-7)
SYNONYMS:	OP-1
RECONSTITUTION:	The lyophilized powder is stable for several weeks at room temperature. For long term storage, reconstitute and aliquot in working volumes. Store reconstituted protein at -20°C.
SOURCE:	CHO cells
PURITY:	> 98% by SDS-PAGE analysis
ACTIVITY:	Induces alkaline phosphatase production in ATDC-5 cells with an ED50 of 0.02 to 0.04 mg/mL.