

Product Specification Sheet

Description: Anti-human VEGF
monoclonal antibody
Product Number: CK55
Aliquot: 1.0 mg
Lot Number: Typical Lot
Storage: -80°C

DESCRIPTION: Monoclonal antibody against Human Vascular Endothelial Growth Factor (VEGF)

MOLECULAR WEIGHT: 38.2 kDa

STORAGE BUFFER: 0.5 x PBS, pH 7.2

STORAGE: Once reconstituted, aliquot and store at -80°C.

IMMUNOGEN: Recombinant Human VEGF

PURIFICATION: IgG1 antibody was purified by Protein G affinity chromatography.

CLASS: IgG1

APPLICATIONS: ELISA 0.2 ng/well.

SPECIES SPECIFICITY: Human.