



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Human Antithrombin
Product Number: CT33
Aliquot: 2.0 mL of 1:1 Slurry
Lot Number: Typical Lot
Storage: 4°C

DESCRIPTION:	Immobilized Human Antithrombin. May be used repeatedly.
MOLECULAR WEIGHT:	58,200
CONCENTRATION:	1.0 mg/mL coupled
SIZE:	1.0 mL
VOLUME:	2.0 mL of a 1:1 slurry
STORAGE:	4°C
FORM:	Frozen liquid
SOURCE:	Human plasma
BUFFER:	0.1M Tris-HCl; 0.15M NaCl; 0.02% Sodium Azide; pH 7.4

For research purposes only.