

Product Specification Sheet

Description: Human LMW
Urokinase, HRP conjugate
Product Number: PA86
Aliquot: 0.1 mg
Lot Number: pa86.typical lot
Storage: -70°C

DESCRIPTION: This Human LMW Urokinase is cleaved from the amino terminal at Lys 135 and retains the proteinase domain. It has been conjugated to HRP.

MOLECULAR WEIGHT: 33 kDa

CONCENTRATION: 2.5 mg/mL

STORAGE BUFFER: 0.05 M Sodium Acetate; 0.1 M NaCl; 1 mM EDTA; pH 5.0

STORAGE: Store at -70°C

SOURCE: Recombinant human urokinase produced in a stable insect cell clone.