

## Human Prostaglandin Synthase-2 (COX-2) Synthetic Peptide

**Product # : PG 97**

**Lot # : 021500YY**

- Peptide:** Distinct C-terminal region of Prostaglandin Synthase-2
- Application:** Will compete with the polyclonal, anti-peptide PGHS-2 antibody in Western blotting, for confirmation of the PGHS-2 specificity of antibody binding.
- Contents:** 100 µg of peptide in 100 µL of 0.1M ammonium bicarbonate buffer, pH 7.5.
- Storage:** Stable at 4°C. For long term storage, aliquot product into separate tubes and freeze at -20°C. Avoid repeated freeze-thaw cycles.