

**Product Specification Sheet**

Description: Anti-Rabbit PAI-1,  
Goat Polyclonal Antibody  
Product Number: PI73  
Aliquot: 1.0 mg  
Lot Number: pi73.604  
Storage: -70°C

<b>DESCRIPTION:</b>	Anti-Rabbit PAI-1 Goat Polyclonal Antibody, IgG Fraction
<b>MOLECULAR WEIGHT:</b>	43 kDa
<b>CONCENTRATION:</b>	8.8 mg/mL
<b>STORAGE BUFFER:</b>	0.05M sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
<b>STORAGE:</b>	-70°C
<b>IMMUNOGEN:</b>	Recombinant full length, wt, non-glycosylated rabbit protein
<b>PURIFICATION:</b>	IgG Fraction
<b>CLASS:</b>	IgG
<b>APPLICATIONS:</b>	EIA, WB
<b>SPECIES SPECIFICITY:</b>	Rabbit