



OXFORD BIOMEDICAL RESEARCH

## Product Specification Sheet

Anti-Protein Kinase C eta ( $\eta$ )

Product Number: PK 45

Aliquot: 100  $\mu$ L

Lot Number 031203DM

Storage: -20°C

<b>DESCRIPTION:</b>	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C- $\eta$ (PKC $\eta$ ).
<b>MOLECULAR WEIGHT:</b>	67.7 kDa Calculated
<b>CONCENTRATION:</b>	Not Determined
<b>PRESENTATION:</b>	Antiserum with 0.05% Sodium Azide
<b>STORAGE:</b>	-20°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
<b>IMMUNOGEN:</b>	KLH conjugated synthetic peptide corresponding to amino acids 666- 683 of PKC $\eta$ . <b>Sequence:</b> <sup>666</sup> N-Q-D-E-F-R-N-F-S-Y-V-S-P-E-L-Q-P <sup>683</sup>
<b>PURIFICATION:</b>	Not Purified
<b>APPLICATIONS:</b>	Western Blotting: 1:1,000 dilution using ECL detection technology.
<b>SPECIES SPECIFICITY:</b>	Human, mouse and rat among other mammalian species.