



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Protein Kinase C theta (θ)

Product Number: PK 55

Aliquot: 100 μ L

Lot Number 020296YY

Storage: -80°C

DESCRIPTION:	This polyclonal antibody is produced in rabbit and is provided as antiserum specific for Protein Kinase C θ (PKC θ).
CONCENTRATION:	Not Determined
PRESENTATION:	Antiserum with 0.05% Sodium Azide
STORAGE:	-80°C. AVOID REPEATED FREEZE THAW CYCLES. Has demonstrated stability at 4 °C for greater than 1 month.
IMMUNOGEN:	KLH conjugated synthetic peptide corresponding to amino acids found on the C-terminus of PKC θ . Sequence: H-S-N-F-S-F-I-N-P-G-M-E-T-L-I-C-S
PURIFICATION:	Not Purified
APPLICATIONS:	Western Blotting: 1:250dilution using alkaline phosphatase detection technology.
SPECIES SPECIFICITY:	Human, mouse and pig among other mammalian species.