



OXFORD BIOMEDICAL RESEARCH

Product Specification Sheet

Anti-Phosphotyrosine IgG

Product Number: PK 99

Aliquot: 250 µg

Lot Number: 020996SS

Storage: -70°C

DESCRIPTION:	This goat polyclonal antibody is specific for phosphotyrosine
MOLECULAR WEIGHT:	The MW of bands on a Western blot is dependent upon the protein being probed.
CONCENTRATION:	0.5 mg/mL as determined using BCA assay with BSA as a reference standards.
STORAGE BUFFER:	PBS
STORAGE:	-70°C. AVOID REPEATED FREEZE THAWS.
IMMUNOGEN:	Phosphotyrosine
CLASS:	IgG
PURIFICATION:	Ammonium Sulfate Fractionation
APPLICATIONS:	Suitable for use in Western blotting, ELISA and immuno-precipitation applications.
SPECIES SPECIFICITY:	Not species specific. May not exhibit utility when using in a goat experimental model.