

Product Specification Sheet

Description: Human PAI-1,
Glycosylated Active Form
Product Number: PL76
Aliquot: 0.5 mg
Lot Number: pl76.915
Storage: -70°C

DESCRIPTION:	Glycosylated wild type human PAI-1 produced in Drosophila S2 cells.
MOLECULAR WEIGHT:	46 kDa
CONCENTRATION:	1.86 mg/ml
STORAGE BUFFER:	0.05M Sodium Phosphate; 0.1M NaCl; 1mM EDTA; pH 6.6
STORAGE:	Store at -70°C
PURIFICATION:	Separated from latent PAI-1 by hydrophobic interaction chromatography.
PURITY:	>95% by SDS-PAGE analysis