



OXFORD BIOMEDICAL RESEARCH

P.O. Box 522, Oxford MI 48371 • USA

USA: 800-692-4633 • Fax: 248-852-4466

www.oxfordbiomed.com

Product specifications

Monoclonal Anti-Rat Cytochrome P450 (2B1, 2B2)

Product # PM 20 Typical Lot

Immunoglobulin

Type: Monoclonal Mouse IgG.

Specificity: Identification of 2B1 and 2B2 by immunoassay or Western blots.

Applications: Immunoassay or Western blots.

Contents: 1mL per vial. Contains sufficient purified monoclonal mouse IgG for at least 10 Western blots.

Working Dilution: A dilution of 1/100 for Western blots will yield a strong band with 10 µg of Phenobarbital induced rat microsomes using Alkaline Phosphatase.

Storage: Stable for several weeks at 4°C, several months at - 20°C, and indefinitely at - 80°C. Do not freeze and thaw repeatedly. Refreezing the product may result in diminishment or loss of activity.