

Product Specification Sheet

Description: Porcine Pancreatic
Elastase, Immobilized
Product Number: SP31
Aliquot: 1.0 mL
Lot Number: sp31.1113
Storage: 4°C

DESCRIPTION:	Immobilized Porcine Pancreatic Elastase. Cleaves PAI-1 to make it non-inhibitory towards uPA and tPA. Coupled to resin and provided as 2 mL of 1:1 slurry.
PURITY:	>95% by SDS-PAGE analysis
SOURCE:	Porcine pancreas
CONCENTRATION:	1.76 mg/mL coupled
STORAGE:	Store at 4°C
BUFFER:	0.1 M HEPES; 0.1 M NaCl; 0.02% Sodium Azide; pH 7.4